# Korrektursocken

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# Corrective socks

## Skewed flatfoot

### Instructions for use of the skewed flatfoot corrective socks

#### How do I put the taping socks on correctly?

Please make sure you use the correct sock for the right and left foot. The 'Taping-Socks' lettering must always be legible on the outer side of the foot.

The taping bands are divided into two function areas: the heel and the ball of the foot.

On the **outer side**, the rear band runs from the base of the Achilles tendon, beneath the outer side of the ankle to the sole of the foot. Its upper edge lightly touches the outer side of the ankle **1**. When pulling the sock on, the band becomes slightly tensioned and thus stabilises the heel in the right position. The front band runs from the sole of the foot at a 45° angle to the ball area at the back of the foot. The front edge sits behind the little toes **2**.

On the **inner side**, the band also runs from the base of the Achilles tendon, beneath the inner side of the ankle towards the sole of the foot. Here, too, its upper edge lightly touches the inner side of the ankle **3**. The tension of the taping band stabilises the heel in its correct position and the longitudinal arch is raised **4**. The front band holds the first metatarsal towards the ground **5**.

#### How do I check the positioning is correct?

The lower leg should be vertically aligned to the foot. The correct positioning of the bands is determined using the 'two finger technique'. For this purpose, the index finger and middle finger are placed around the instep **6**. If the upper edge of the front band is in line with the middle finger, it is correctly positioned **7**.

Outer side

Inner side